

## Mineral Insulated Dual Thermocouple with Standard Duplex Type K Plug



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

### Need Your Temperature Sensor calibrated?

Traceable to ISO17025 (UKAS) calibration standards, Labfacility offer a choice of a 3 or 5 point traceable calibration using our in-house calibration facility.

**Please be aware that if you add a calibration option to your order you will be contacted regarding delivery information.**

### Specifications

## Specifications

|                              |  |
|------------------------------|--|
| <b>General Description</b>   | 310 stainless steel sheath. Sheath can be bent/formed to suit many applications and processes. |
| <b>Sensor Type</b>           | Mineral insulated Type 'K' with Duplex Thermocouple Plug                                       |
| <b>Thermocouple Type</b>     | K IEC  |
| <b>Thermocouple Junction</b> | Insulated hot junction   |
| <b>Termination Type</b>      | Standard duplex plug termination (200°C)   |
| <b>Max. Temperature</b>      | 1100°C   |
| <b>Min. Temperature</b>      | -40°C  |
| <b>Standards Met</b>         | Conforms to IEC 584 specification  |
| <b>Connector Type</b>        | Plug   |
| <b>Connector Size</b>        | Standard   |
| <b>Length</b>                | 150, 250 or 500mm  |