

## Handheld Needle Probe T-shaped - Asphalt or Food Processing Type K - 1M length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

T-shaped handle reduced tip penetration probe.

Reduced sensing tip for faster response.

Response time less than 17 seconds.

Probe temperature range -75 to 250°C

Accuracy to BS EN 60584-1:2013

Tip diameter 9.5 x 1000mm or 1400mm length

Biomaster antimicrobial technology is incorporated during the manufacturing process, ensuring lifelong protection against the risk of bacterial cross contamination. As well as being harmful to health, outbreaks of food poisoning can damage brands, and the reputation of stores, food producers and restaurants.

The risk of food-borne contamination is widespread, existing in restaurants, bars, hotels, supermarkets, hospitals, schools, aircraft, cruise ships and your home. Anywhere where food is made, prepared or eaten.

### Specifications

## Specifications

<b>Product Code</b>	XE-3798-001
<b>General Description</b>	Probe response time: less than 17 seconds
<b>Sensor Type</b>	T-shaped handle reduced tip penetration probe
<b>Thermocouple Type</b>	K IEC
<b>Probe Diameter</b>	9.5mm
<b>Probe Length</b>	1000mm
<b>Max. Temperature</b>	250°C
<b>Min. Temperature</b>	-75°C
<b>Standards Met</b>	BS EN 60584-1:2013