



## Handheld Surface Probe - Type K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Stainless steel sheath with a moulded handle and 2 metres of coiled extension cable terminated in a miniature thermocouple plug.

## Specifications

<b>Product Code</b>	XE-3840-001
<b>General Description</b>	Handheld Surface Probe
<b>Thermocouple Type</b>	K
<b>Calibration</b>	IEC
<b>Probe Diameter</b>	15mm
<b>Probe Length</b>	170mm
<b>Sheath Material</b>	Stainless steel sheath
<b>Cable Length</b>	2m
<b>Cable Termination</b>	Fitted Plug
<b>Max. Temperature</b>	400°C