

## Handheld Moving Air Probe - Type K,T



[Click here to view product](#)

<b>Product Code</b>	XE-3845
<b>General Description</b>	Handheld Moving Air Probe, Ideal for the measurement of air and gases
<b>Sensor Type</b>	Funnelled shield allows for maximum airflow across measuring junction
<b>Thermocouple Type</b>	K or T
<b>Calibration</b>	IEC
<b>Probe Diameter</b>	30mm
<b>Probe Length</b>	120mm
<b>Cable Length</b>	2m
<b>Cable Termination</b>	Fitted Plug