

Handheld Moving Air Probe - Type K IEC 30mm x 120mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Stainless steel sheath with a moulded handle and 2 metres of coiled extension cable terminated in a miniature thermocouple plug. Thermocouple air probes are ideal for measuring the temperature of still or moving air and gases.

Specifications

Specifications

Product Code	XE-3845-001
General Description	Handheld Moving Air Probe, Ideal for the measurement of air and gases
Sensor Type	Funnelled shield allows for maximum airflow across measuring junction
Thermocouple Type	K IEC
Probe Diameter	30mm
Probe Length	120mm
Cable Length	2m
Cable Termination	Fitted plug
Series	Handheld Probe