

## Handheld Moving Air Probe Type T IEC 30mm Dia x 120mm Length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Stainless steel sheath with a moulded handle and 2 metres of coiled extension cable terminated in a miniature thermocouple plug. Thermocouple air probes are ideal for measuring the temperature of still or moving air and gases.

### Specifications

## Specifications

<b>Product Code</b>	XE-3846-001
<b>General Description</b>	Handheld Moving Air Probe, Ideal for the measurement of air and gases. Funnelled shield allows for maximum airflow across measuring junction
<b>Thermocouple Type</b>	T IEC
<b>Probe Diameter</b>	30mm
<b>Probe Length</b>	120mm
<b>Cable Length</b>	2m
<b>Cable Termination</b>	Fitted plug