


## Pipe Clamp Probe - type K



[Click here to view product](#)

<b>Product Code</b>	XE-3847-001
<b>General Description</b>	Pipe Clamp Probe
<b>Sensor Type</b>	Clamps onto pipes from 6 mm to 30 mm
<b>Thermocouple Type</b>	K
<b>Cable Length</b>	1m
<b>Cable Type</b>	PVC
<b>Cable Termination</b>	Miniature Plug
<b>Max. Temperature</b>	-10°C
<b>Min. Temperature</b>	+100°C
<b>Accuracy</b>	Response time less than 2 seconds