

## Pipe Clamp Probe - type K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

This pipe clamp probe is suitable for measuring the surface temperature of pipes in heating and ventilating systems, refrigeration and many other pipe applications.

## Specifications

<b>Product Code</b>	XE-3847-001
<b>General Description</b>	Pipe Clamp Probe
<b>Sensor Type</b>	Clamps onto pipes from 6 mm to 30 mm
<b>Thermocouple Type</b>	K
<b>Cable Length</b>	1m
<b>Cable Type</b>	PVC
<b>Cable Termination</b>	Miniature Plug
<b>Max. Temperature</b>	-10°C
<b>Min. Temperature</b>	+100°C
<b>Accuracy</b>	Response time less than 2 seconds