

Pipe Clamp Probe - type K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

This pipe clamp probe is suitable for measuring the surface temperature of pipes in heating and ventilating systems, refrigeration and many other pipe applications.

Specifications

Product Code	XE-3847-001
General Description	Pipe Clamp Probe
Sensor Type	Clamps onto pipes from 6 mm to 30 mm
Thermocouple Type	K
Cable Length	1m
Cable Type	PVC
Cable Termination	Miniature Plug
Max. Temperature	-10°C
Min. Temperature	+100°C
Accuracy	Response time less than 2 seconds