

## Handheld Stainless Steel Penetration Probe - PRT (black) 3.3mm x 120mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

<p>316 Stainless Steel Penetration Probe with a pointed tip for liquid and semi-solid temperature measurement .  
 Comprising of a handle, 2 metre coiled cable and a mini plug.</p>

<p> </p>  
 <table style="height: 36px; border-color: #f5f5f5;" border="1" width="850">

<tbody>  
 <tr>

<th style="width: 417px; background-color: #cccccc; text-align: center; vertical-align: middle;">  
 <h2><strong>Need Your Temperature Sensor calibrated?</strong></h2>

</th>

<td style="width: 417px; background-color: #f5f5f5; text-align: center; vertical-align: middle;">

<p>Traceable to ISO17025 (UKAS) calibration standards, Labfacility offer a choice of a 3 or 5 point traceable calibration using our in-house calibration facility, you can select any temperature points between -10°C & +200°C.</p>

<p><strong>Please be aware that if you add a calibration option to your order you will be contacted regarding delivery information.</strong></p>

</td>

</tr>

</tbody>

</table>

<p><span style="text-decoration: underline;"><strong>Specifications</strong></span></p>

<p> </p>

## Specifications

<b>Product Code</b>	XE-3890-001
<b>General Description</b>	Penetration Probe with a pointed tip
<b>Sensor Type</b>	PRT
<b>Probe Diameter</b>	3.3mm
<b>Probe Length</b>	120mm
<b>Cable Length</b>	2m
<b>Max. Temperature</b>	+400°C