


GM520 Pressure Manometer



[Click here to view product](#)

Product Code	XE-4011-001
Model	GM520 Pressure Manometer
Communications Protocol	USB communication
Accuracy	+/- 0.3% FSO (25Å°C)
Repeatability	+/- 0.2% (Maximum +/- 0.5% FSO)
Response	Typical 0.5 s
Linearity/hysteresis	+/- 0.29% FSO
Maximum Pressure	+/- 150kPa
Storage Temperature Range	-10 to 60Å°C
Ambient Operating Temperature	0 to 50Å°C

Measuring range	+/- 35kPa
-----------------	-----------

