

Adjustable Racing Tyre Temperature Probe



[Click here to view product](#)

Product Code	XE-3892-001
General Description	This very fast response racing tyre probe has a manual depth adjustment from 1 mm to 15 mm
Sensor Type	Can be used with any type K digital thermometer with a mini input socket.
Thermocouple Type	K
Probe Diameter	1.1mm
Probe Length	manual depth adjustment from 1 mm to 15 mm
Cable Length	1m
Cable Termination	Miniature Plug
Max. Temperature	+250°C

Min. Temperature	-50°C
Features #1	Junction Type: Grounded

