

## Handheld Between Pack / Flat Food Probe - Type K,T



[Click here to view product](#)

<b>Product Code</b>	XE-3848
<b>Sensor Type</b>	Hand Held Between Pack/Flat Food Probe
<b>Thermocouple Type</b>	K or T
<b>Calibration</b>	IEC
<b>Thermocouple Junction</b>	Single element, grounded junction offering fast response.
<b>Probe Diameter</b>	4mm
<b>Probe Length</b>	140mm approximately below the handle
<b>Cable Length</b>	2 metre (extended)
<b>Cable Type</b>	Polyurethane
<b>Cable Termination</b>	moulded on miniature plug.

<b>Max. Temperature</b>	+260°C
<b>Min. Temperature</b>	-50°C
<b>Standards Met</b>	IEC 584

