


Handheld Between Pack / Flat Food Probe - Type K,T



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

Sensor Type	Hand Held Between Pack/Flat Food Probe
Thermocouple Type	K or T
Calibration	IEC
Thermocouple Junction	Single element, grounded junction offering fast response.
Probe Diameter	4mm
Probe Length	140mm approximately below the handle
Cable Length	2 metre (extended)
Cable Type	Polyurethane
Cable Termination	moulded on miniature plug.
Max. Temperature	+260°C
Min. Temperature	-50°C
Standards Met	IEC 584