

Digital Catering Thermometer Type T



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Catering Thermometer

High accuracy chef thermometer

Thermocouple type T

HACCP compliant,

-100 to + 280°C

Hold function,

IP67 waterproof thermometer

Dual LCD display

Easy to use: two button operation

Dual display

Auto switch off capability.

Features:

Self-calibrating

Low battery indication

Supplied with 1 9V PP3 battery

Measurement Range -100 to +280°C

High Accuracy +/-0.1% of Reading

+/-0.2°C, 0.1°C Resolution

The CA2005 offers a hold function and when this button is pressed, the current temperature will be held and transferred to the top right of the display. The main display will continue to show the temperature in real-time.

Specifications

Specifications

Product Code	XE-4090-001
General Description	High accuracy chef thermometer
Type	IP67 waterproof thermometer

Accuracy	+/-0.1% of Reading
Resolution	+/-0.2°C, 0.1°C Resolution
Input / Ranges	-100 to + 280°C
Power Supply	Supplied with 1 9V PP3 battery
Sensor Type	T