



## **PTFE Single Shot Thermocouple Cable Type K**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**PTFE Single Shot Thermocouple Cable Type K**

Unique single shot cable construction Resistant to chemicals, oils and fluids. Suitable for aerospace applications and is widely used in autoclaves.

**Specifications****Specifications**

<b>General Description</b>	Unique single shot cable construction Resistant to chemicals, oils and fluids. Suitable for aerospace applications and is widely used in autoclaves.
<b>Cable Length</b>	10, 25, 50, 100M
<b>Cable Type</b>	PTFE
<b>Core / Strands</b>	1/0.376mm
<b>Colour</b>	Green
<b>Max. Temperature</b>	+250°C
<b>Min. Temperature</b>	-75°C
<b>Standards Met</b>	IEC-584-3
<b>Thermocouple Type</b>	K