

PTFE Single Shot Thermocouple Cable Type T



50M



10m



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

PTFE Single Shot Thermocouple Cable Type T

Unique single shot cable construction Resistant to chemicals, oils and fluids. Suitable for aerospace applications and is widely used in autoclaves.

Specifications**Specifications**

General Description	Unique single shot cable construction Resistant to chemicals, oils and fluids. Suitable for aerospace applications and is widely used in autoclaves.
Cable Length	10, 25, 50, 100M
Cable Type	PTFE
Core / Strands	1/0.376mm
Colour	Brown
Max. Temperature	+250°C
Min. Temperature	-75°C
Standards Met	IEC-584-3, Class 1
Thermocouple Type	T