

PFA Insulated Flat Pair IEC Exposed Junction Thermocouple with Miniature Plug - Type K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Good mechanical strength and flexibility, resistant to oils, acids and other adverse fluids.

Need Your Temperature Sensor calibrated?

Traceable to ISO17025 (UKAS) calibration standards, Labfacility offer a choice of a 3 or 5 point traceable calibration using our in-house calibration facility, you can select any temperature points between -10°C & +200°C.

Please be aware that if you add a calibration option to your order you will be contacted regarding delivery information

Specifications

Sensor Type	Fast response, welded exposed junction
Thermocouple Type	Thermocouple type K
Cable Connection	Terminated with a miniature flat pin plug
Cable Length	300mm, 1, 2, 3 or 5 metre lengths available
Core / Strands	1/0.3mm or 1/0.5mm
Termination Type	Miniature Plug
Standards Met	(cores & jacket) to IEC-584-3
Insulation Material	Insulation rating: -75°C to 260°C