

## Terminal Blocks - Type K IEC



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Terminal blocks come in many different types and sizes to suit your terminal head.

Check carefully the dimensions of the fixing holes on your terminal head before choosing the correct block.

**Way**

**Order Code**

**Colour**

**Diameter**

**Fixings**

**Fits Head**

2

XE-3825-001

Green

42mm

33mm

KNE, KNP, KPP, B, SCH4

4

XE-3826-001

Green

42mm

33mm

KNE, KNP, KPP, B, SCH4

## Specifications

<b>Material</b>	Plastic
<b>Fixing (Centres)</b>	33mm
<b>Colour</b>	Green
<b>Dimensions</b>	42mm Diameter
<b>To Fit Head</b>	KNE, KNP, KPP, B, SCH4
<b>Sizes</b>	2 or 4 Way
<b>Thermocouple Type</b>	K