

## KNP Hinged Lid Terminal Head



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Many alternative types of terminal heads are available to meet the requirements of various applications.

Variations exist in size, material, accommodation, resistance to media, resistance to fire or even explosion and in other parameters.

## Specifications

<b>Head Type</b>	KNP Hinged Lid Terminal Head
<b>Material</b>	Aluminium Silver Painted
<b>Cable Entry</b>	M20 x 1.5mm
<b>Fixing (Centres)</b>	33mm & 40mm
<b>Max. Temperature</b>	170°C
<b>IP Rating</b>	IP67 rated
<b>Process Entry</b>	1/2" BSPP and 3/4" BSPP
<b>Cap Retaining Type</b>	Hinged Cap