

Legionella Thermometer with Integral Water Temp Timer



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

MM2008 Single input Thermometer with an integral 1 minute & 2 minute timer.

This Thermometer is primarily developed for use in Legionella risk management and offers reassurance that the correct temperature is taken when monitoring hot and cold water temperatures.

Features

Single Input Thermocouple

Integral Timer for 1 and 2 minute countdown

Measurement Range: -200 to +1,372°C (probe type dependant)

Accuracy: - $\pm 0.1\%$ of Reading $\pm 0.2^\circ\text{C}$

C/F Selectable

IP67 Waterproof Casing

Low Battery Indication

Selectable Thermocouple Types: K & T

Self-Calibrating

Rubber Boot included for extra protection against damage

Includes 1 x PP3 Battery

Requires: A Thermocouple Probe/Sensor

Specifications

Specifications

Product Code	XE-4091-001
General Description	Measurement Range: -200 to +1,372°C (probe type dependant)
Type	Single input Thermometer with an integral 1 minute & 2 minute timer
Accuracy	$\pm 0.1\%$ of Reading $\pm 0.2^\circ\text{C}$
Readings	C/F Selectable
Power Consumption	Low Battery Indication
Power Supply	Includes 1 x PP3 Battery



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Standards Met	IP67 Waterproof Casing
Thermocouple Type	K & T
Thermometer	Self-Calibrating
Termination Type	Requires a Thermocouple Probe/Sensor
Accessories	Rubber Boot included for extra protection against damage