


Type J Glassfibre Insulated Flat Pair thermocouple Cable / Wire (BS)

Type J, 1/0.508mm
Fiberglass Flat pair



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Type J Glassfibre Insulated Flat Pair thermocouple Cable / Wire (BS)

Image is for colour coding illustration only

Price is per metre - please add required quantity for basket. Only 1,227 metres left in stock

Glassfibre insulated flat pair construction, conductors laid flat, Glassfibre insulated with Glassfibre overall, silicone varnished throughout. Good temperature resistance but not suitable where fluids are present.

Specifications
Specifications

Product Code	WJ-033-D
General Description	Good temperature range but will not prevent ingress of fluids. Fairly flexible but does not provide good mechanical protection.
Cable Length	Priced per metre
Cable Type	Flat pair cable
Core / Strands	1/0.508mm
Outer Colour	Black
+ Colour	Yellow
- Colour	Blue
Insulation Material	Glassfibre insulated
Max. Temperature	+350°C (short periods up to 400°C)
Min. Temperature	-60°C
+ Conductor	Iron
- Conductor	Constantan
Standards Met	BS



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Thermocouple Type	J
Tolerance	Class 2