


## Type J Glassfibre Insulated Flat Pair thermocouple Cable / Wire (BS)

Type J, 1/0.508mm  
Fiberglass Flat pair



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Type J Glassfibre Insulated Flat Pair thermocouple Cable / Wire (BS)**

Image is for colour coding illustration only

**Price is per metre - please add required quantity for basket. Only 1,227 metres left in stock**

Glassfibre insulated flat pair construction, conductors laid flat, Glassfibre insulated with Glassfibre overall, silicone varnished throughout. Good temperature resistance but not suitable where fluids are present.

**Specifications**
**Specifications**

<b>Product Code</b>	WJ-033-D
<b>General Description</b>	Good temperature range but will not prevent ingress of fluids. Fairly flexible but does not provide good mechanical protection.
<b>Cable Length</b>	Priced per metre
<b>Cable Type</b>	Flat pair cable
<b>Core / Strands</b>	1/0.508mm
<b>Outer Colour</b>	Black
<b>+ Colour</b>	Yellow
<b>- Colour</b>	Blue
<b>Insulation Material</b>	Glassfibre insulated
<b>Max. Temperature</b>	+350°C (short periods up to 400°C)
<b>Min. Temperature</b>	-60°C
<b>+ Conductor</b>	Iron
<b>- Conductor</b>	Constantan
<b>Standards Met</b>	BS



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Thermocouple Type	J
Tolerance	Class 2