


Type J Glassfibre Insulated Flat Pair Cable / Wire with Stainless Steel Overbraid (BS)

Type J, 14/0.2mm
SSOB PTFE



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Type J Glassfibre Insulated Flat Pair Cable / Wire with Stainless Steel Overbraid (BS)

Image is for colour coding illustration only

Price is per metre - please add required quantity for basket. Only 336 metres left in stock

Glassfibre insulated flat pair construction, conductors laid flat, Glassfibre insulated with Glassfibre jacket, silicone varnished throughout + stainless steel overbraid for added protection and mechanical strength. Good temperature resistance but not suitable where fluids are present.

Specifications
Specifications

Product Code	WJ-247-D
General Description	Good temperature range but will not prevent ingress of fluids. Fairly flexible but does not provide good mechanical protection.
Cable Length	Priced per metre
Cable Type	SSOB PTFE
Core / Strands	14/0.2mm
Outer Colour	Black
+ Colour	Yellow
- Colour	Blue
Insulation Material	Glassfibre insulated extension cable with stainless steel overbraiding
Max. Temperature	+350°C (short periods up to 400°C)
Min. Temperature	-60°C
+ Conductor	Iron
- Conductor	Constantan



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Standards Met	BS
Thermocouple Type	J
Tolerance	Class 2