


# IEC Thermocouple Brake Pad Temperature Sensors Type K, J



[Click here to view product](#)

<b>Product Code</b>	XE-3893
<b>General Description</b>	Designed for monitoring the temperature of brake pads
<b>Sensor Type</b>	Grounded Junction
<b>Thermocouple Type</b>	K or J
<b>Probe Diameter</b>	3.2mm, 4.0mm or 6.0mm
<b>Probe Length</b>	Tip Length: 4mm
<b>Cable Length</b>	1m
<b>Cable Type</b>	Fibre glass with stainless steel Overbraid 7/0.2mm
<b>Termination Type</b>	Bare tails
<b>Max. Temperature</b>	+400°C

Features #1	Potting sealant: Ceramic paste (Fortafix)
-------------	---

