

Pt100 Duplex 2 x 3 wire class A Resistance Thermometer, KNE Head



[Click here to view product](#)

Product Code	XE-3655
Sensor Type	2 x Pt100 (100 Ohms @ 0°C) to IEC 751, Class A
Probe Diameter	6mm
Probe Length	100mm, 150mm, 200mm or 250mm
Sheath Material	316 stainless steel
Head Type	KNE
Termination Type	6 way connection block for KNE head
Max. Temperature	+250°C (connection head @ 170°C)
Min. Temperature	-50°C
IP Rating	IP67

Accuracy	Class A
Number of Wires	2 x 3 Wire

