

Industrial Style Mineral Insulated Thermocouple with KNE Terminal Head - Type K Inconel (6mm x 400mm)



[Click here to view product](#)

Product Code	XE-3248-001
General Description	Highly flexible, sheath can be bent/formed to suit many applications and processes
Sensor Type	Mineral insulated Thermocouple fitted with KNE terminal head
Thermocouple Type	Type K
Thermocouple Junction	Insulated hot junction
Probe Diameter	6.0mm diameter
Probe Length	400mm Length
Sheath Material	Inconel
Head Type	KNE, Cast aluminium finish, screw type lid with chain
Gland	M16 x 1.5mm Plated brass gland fitted

Max. Temperature	+1100°C
Min. Temperature	-40°C
IP Rating	IP67

