

Digital Body Thermometer



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

The digital thermometer can quickly, accurately and conveniently measure the temperature of the human body. The shell is made of non-toxic plastic and is less fragile than the traditional mercury thermometer, which is easy to break. Mercury pollutes the environment and is not as easy to read.

Specifications

Specifications

Product Code	XE-4095-001
General Description	Digital Thermometer
Type	Flexible probe head
Accuracy	+/- 0.1°C
Resolution	0.1°C
Dimensions	124mm x 18.5mm x 10mm
Temperature Range	+32°C to 42°C
Weight	10g
Important Information	Measurement time: about 1 minute