




## Digital Body Thermometer



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

The digital thermometer can quickly, accurately and conveniently measure the temperature of the human body.

The shell is made of non-toxic plastic and is less fragile than the traditional mercury thermometer, which is easy to break. Mercury pollutes the environment and is not as easy to read.

### Specifications

## Specifications

<b>Product Code</b>	XE-4095-001
<b>General Description</b>	Digital Thermometer
<b>Type</b>	Flexible probe head
<b>Accuracy</b>	+/- 0.1°C
<b>Resolution</b>	0.1°C
<b>Dimensions</b>	124mm x 18.5mm x 10mm
<b>Temperature Range</b>	+32°C to 42°C
<b>Weight</b>	10g
<b>Important Information</b>	Measurement time: about 1 minute