

0.25mm Diameter Fast Response Mineral Insulated Thermocouple with ANSI Miniature Plug - Type K



IMAGE
COMING
SOON

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Please use our [custom sensor forms](#) to submit your requirements for various lengths of MI, as all terminated with an IEC miniature plug (220°C simple connection).

Need Your Temperature Sensor calibrated?

Traceable to ISO17025 (UKAS) calibration standards, Labfacility offer a choice of a 3 or 5 point traceable calibration using our in-house calibration facility.

Please be aware that if you add a calibration option to your order you will be contacted regarding delivery information.

Specifications

Specifications

General Description	Highly flexible, sheath can be bent/formed to suit many applications and processes
Sensor Type	Type K ANSI MB STYLE probes 0.25mm
Thermocouple Junction	Single element, junction insulated or grounded
Probe Length	100mm / 150mm / 250mm
Sheath Material	Stainless steel 316 sheath
Termination Type	Fitted Miniature Plug Type ANSI
Range	-40 to +260°C
Important Information	A 1%+ Length Tolerance applies