

## Self Adhesive Patch PT100 Sensor with PFA insulated outer jacket - Type PRT



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

For more options including thermistor type, RTD detector type, varying detector grades and custom lengths, please contact our sales team using the product enquiry tab.

The RTF4/NAH style is a surface sensor encapsulated within a silicone rubber patch.

Used for surface applications

Class B tolerance

PFA insulated lead

4-wire tail termination to IEC60751

1, 3 or 5m PFA lead options available

### Specifications

## Specifications

<b>Sensor Type</b>	Silicone rubber patch with self adhesive foil backing
<b>Element</b>	Pt100
<b>Patch dimensions</b>	L40 x W15mm x D5mm
<b>Cable Length</b>	1, 3 or 5m
<b>Cable Type</b>	PFA insulated outer jacket
<b>Core / Strands</b>	7/0.1mm
<b>Cable Termination</b>	Bare Tails
<b>Max. Temperature</b>	+150°C
<b>Min. Temperature</b>	-50°C
<b>Accuracy</b>	Class B
<b>Number of Wires</b>	4