

0.5mm Diameter with threaded pot & tails ? Insulated or Grounded



[Click here to view product](#)

| | |
|-----------------------------|--|
| Product Code | XE-0154 |
| General Description | Type K with threaded pot & tails in 310 stainless steel |
| Sensor Type | Mineral Insulated |
| Thermocouple Type | K |
| Element/hot junction | Single element, junction insulated or grounded from sheath |
| Probe Diameter | 0.5mm |
| Probe Length | 50, 100, 150, 250 mm |
| Sheath Material | 310 stainless steel |
| Pot Type | M8 x 1.0mm fine pitch threaded pot seal (200°C) |
| Cable Length | 100mm |

| | |
|-------------------------|--------------------------------|
| Cable Type | PFA insulated twin twisted IEC |
| Core / Strands | 7/0.2mm |
| Termination Type | Bare Tails |
| Max. Temperature | +650°C |
| Min. Temperature | -40°C |
| Standards Met | In accordance with IEC 584 |