


## Digital Household Thermometer Hygrometer with rear magnet



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

TypeT-1/0.2

## Specifications

<b>Product Code</b>	XE-4352-001
<b>Model</b>	LDTH-8
<b>General Description</b>	Digital Household Thermometer Hygrometer with rear magnet
<b>Type</b>	Timer-Thermometer 2 in 1
<b>Accuracy</b>	+/- 3%
<b>Readings</b>	Timer range: 99 Min 59 Sec
<b>Resolution</b>	0.1°C(0.2°F)
<b>High Alarm Level</b>	High temperature alarm function
<b>Battery</b>	CR2430
<b>Cable Length</b>	Wire: 120cm length
<b>Length</b>	10m
<b>Colour</b>	White
<b>Dimensions</b>	118mm x 60mm x 7.5mm
<b>Functions</b>	°C/°F
<b>Function #1</b>	Continuous: 33°F- 140°F
<b>Display Type</b>	LCD
<b>Operating Humidity</b>	25%-98% RH
<b>Temperature Range</b>	-50-70°C (-58-158°F)
<b>Probe Material</b>	304 stainless steel
<b>Weight</b>	86g
<b>Sensor Type</b>	Sensirion SHT30



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

<b>Standards Met</b>	CE, FCC, ROHS
<b>Construction</b>	Eco-friendly ABS plastic + Stainless steel front shell
<b>Display</b>	FSLCD
<b>Important Information</b>	Support stand or magnetic mounting options