

Digital folding probe thermometer LDT802



[Click here to view product](#)

Product Code	XE-4351-001
Model	LDT802
General Description	Digital folding probe thermometer
Accuracy	+/-1Â°C (-2Â°F) at -20Â°C to150Â°C (-4Â°F to 392Â°F)
Readings	Min max temperature record
Resolution	0.1Â°C/0.2Â°F
Response	Super fast response time of 3-4 seconds
Battery	3V button battery (included)
Display Type	Bright backlit
Automatic Shutdown	Power saving: 10 minutes auto-off

Temperature Range	-50°C to 300°C (-58°F to 572°F)
Probe Diameter	Foldable probe Diameter: 3.5mm
Probe Material	304 stainless steel
Standards Met	IP67 Waterproof design
Construction	Eco-friendly ABS plastic
Important Information	Inside magnet on backside