

## JIS Type K Exposed Welded Tip Thermocouple 1/0.3mm PFA Twin Twisted Cable with Miniature Plug 2m



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Mini plug termination.

Welded, exposed junction for fast response times. Durable construction with PFA insulation for resistance to oils, acid, fluids and good mechanical strength and flexibility.

## Specifications

<b>Product Code</b>	XE-0520-001
<b>General Description</b>	Welded, exposed junction for fast response times. Durable construction with PFA insulation for resistance to oils, acid, fluids and good mechanical strength and flexibility.
<b>Manufacturing Code</b>	Z2-K-2.0-1/0.3-MP-JIS
<b>Sensor Type</b>	Exposed Junction
<b>Thermocouple Type</b>	JIS Type K
<b>Cable Length</b>	2m (6.56ft)
<b>Cable Type</b>	PFA Twin Twisted
<b>Core / Strands</b>	1/0.3mm
<b>Cable Termination</b>	Mini Plug
<b>Max. Temperature</b>	+250°C
<b>Min. Temperature</b>	-75°C