



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

JIS Type K Exposed Welded Tip Thermocouple 1/0.3mm PFA Twin Twisted Cable with Miniature Plug 5m



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Mini plug termination.

Welded, exposed junction for fast response times. Durable construction with PFA insulation for resistance to oils, acid, fluids and good mechanical strength and flexibility.

Specifications

Product Code	XE-0521-001
General Description	Welded, exposed junction for fast response times. Durable construction with PFA insulation for resistance to oils, acid, fluids and good mechanical strength and flexibility.
Manufacturing Code	Z2-K-5.0-1/0.3-MP-JIS
Sensor Type	Exposed Junction
Thermocouple Type	JIS Type K
Cable Length	5m (16.4ft)
Cable Type	PFA Twin Twisted
Core / Strands	1/0.3mm
Cable Termination	Mini Plug
Max. Temperature	+250°C
Min. Temperature	-75°C