

JIS Type K Exposed Welded Tip Thermocouple 1/0.2mm PFA Twin Twisted Cable with Bare Tails 1m



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Bare Tails Termination.

Welded, exposed junction for fast response times. Durable construction with PFA insulation for resistance to oils, acid, fluids and good mechanical strength and flexibility.

Specifications

Product Code	XE-0511-001
General Description	Welded, exposed junction for fast response times. Durable construction with PFA insulation for resistance to oils, acid, fluids and good mechanical strength and flexibility.
Manufacturing Code	Z2-K-1.0-1/0.2-T-JIS
Sensor Type	Exposed Junction
Thermocouple Type	JIS Type K
Cable Length	1m (3.3ft)
Cable Type	PFA Twin Twisted
Core / Strands	1/0.2mm
Cable Termination	Bare Tails
Max. Temperature	+250°C
Min. Temperature	-75°C