



IEC Type K 4.5kg Pull Power (Horseshoe) Magnet Thermocouple, PFA Insulated Cable with Stainless Steel Over-Braid Terminating in Miniature or Standard Plug



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

General Description	4.5kg Pull Power (Horseshoe) Magnet Thermocouple
Sensor Type	Magnet Thermocouple
Thermocouple Type	Type K
Magnet	30 x 20 x 20mm (L x H x W). Max Pull: 4.5kg
Cable Length	1.0M or 2.0M
Cable Type	PFA insulated lead with stainless steel overbraid
Core / Strands	7/0.2mm
Termination Type	Miniature or Standard Plug
Max. Temperature	250°C
Min. Temperature	-50°C
Standards Met	IEC