


## JIS Type K Magnetic Strip Thermocouple, PFA Insulated with Miniature Plug



[Click here to view product](#)

<b>Product Code</b>	XE-0553
<b>General Description</b>	Ideal for ferrous metal surface measurement
<b>Thermocouple Junction</b>	Grounded measuring junction
<b>Magnet</b>	50 x 25 x 1mm magnetic strip
<b>Cable Length</b>	1 or 2 metre
<b>Cable Type</b>	PFA insulated wires, twin twisted
<b>Termination Type</b>	Miniature flat pin plug
<b>Max. Temperature</b>	100°C
<b>Min. Temperature</b>	-50°C
<b>Standards Met</b>	JIS (Japanese Industrial Standards)

<b>Colour</b>	Positive: Red / Negative: White
---------------	---------------------------------

