




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

JIS Type K 0.7kg Pull Button Magnet Thermocouple, PTFE Insulated with Miniature Plug - 0.5m



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-0550-001
General Description	Ideal for ferrous metal surface measurement
Sensor Type	K JIS
Thermocouple Junction	Grounded at tip
Magnet	1/2" diameter button type, 0.7kg Pull
Cable Length	0.5m
Cable Type	PFA insulated twin-twist wires
Cable Termination	Miniature Fitted Plug
Max. Temperature	+250°C (plug 200°C)
Standards Met	JIS (Japanese Industrial Standards)
Colour	Positive: Red / Negative: White