




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

ANSI Fine Wire Disc Thermocouple, PFA with Miniature Plug - Type K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

General Description	Exposed junction thermocouple with Copper Disc brazed at tip
Sensor Type	Copper disc, 6.35mm diameter x 0.3mm thick
Thermocouple Type	Type K
Thermocouple Junction	Grounded measuring junction
Cable Length	1m, 2m, 5m
Cable Type	1/0.2mm PFA Teflon® insulated twisted wire
Termination Type	Miniature Plug or Bare tails
Max. Temperature	200°C
Min. Temperature	0°C
Standards Met	ANSI MC96.1
Colour	Positive: Yellow / Negative: Red
Accessories	Supplied with high temperature adhesive kapton strips (13x18mm x 5) for quick and easy attachment to surface/process