



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

ANSI Bolt Thermocouple, Glassfibre Stainless Steel Over Braided Cable - Type K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Bolt sensors are generally screwed into the surface which normally increases the speed of response and sensitivity of the measurement. Suitable for use up to +250°C.

Specifications

Specifications

Thermocouple Type	K
Thermocouple Junction	Bolts have grounded junctions
Cable Length	2M
Cable Type	Glassfibre stainless steel over braided cable
Cable Termination	Bare Tails
Max. Temperature	+250°C
Standards Met	ANSI MC96.1
Bolt thread/pitch	M8 x 1mm / M10 x 1mm / M12 x 1mm
Colour	Positive: Red / Negative: Yellow
Length	Bolt length: 13mm or 25mm