

**Bolt Thermocouple, Glassfibre Stainless Steel  
Over Braided Cable - JIS Type K- M8x1mm thread  
pitch - 25mm bolt size**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Bolt sensors are generally screwed into the surface which normally increases the speed of response and sensitivity of the measurement. Suitable for use up to +250°C.

### Specifications

## Specifications

<b>Product Code</b>	XE-0591-001
<b>Thermocouple Type</b>	K JIS
<b>Thermocouple Junction</b>	Bolts have grounded junctions
<b>Cable Length</b>	2m
<b>Cable Type</b>	Glassfibre stainless steel over braided cable
<b>Cable Termination</b>	Bare Tails
<b>Max. Temperature</b>	250°C
<b>Standards Met</b>	JIS (Japanese Industrial Standards)
<b>Bolt thread/pitch</b>	M8 x 1mm
<b>Colour</b>	Positive: Red / Negative: White
<b>Length</b>	25mm