



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Bolt Thermocouple, Glassfibre Stainless Steel Over Braided Cable - JIS Type K- M12x1mm thread pitch - 25mm bolt size



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Bolt sensors are generally screwed into the surface which normally increases the speed of response and sensitivity of the measurement. Suitable for use up to +250°C.

Specifications

Specifications

Product Code	XE-0595-001
Thermocouple Type	K JIS
Thermocouple Junction	Bolts have grounded junctions
Cable Length	2m
Cable Type	Glassfibre stainless steel over braided cable
Cable Termination	Bare Tails
Max. Temperature	250°C
Standards Met	JIS (Japanese Industrial Standards)
Bolt thread/pitch	M12 x 1mm
Colour	Positive: Red / Negative: White
Length	25mm