



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

ANSI Type K 1.9kg Pull Button Magnet Thermocouple, PFA Insulated Cable with Stainless Steel Over-Braid Terminating in Miniature or Standard Plug



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

| | |
|----------------------------|--|
| General Description | 1.9kg Pull Button Magnet Thermocouple |
| Thermocouple Type | K |
| Magnet | 19.1mm Diameter / 12.7mm height |
| Cable Length | 1.0m or 2.0m |
| Cable Type | Glassfibre insulated lead with stainless steel overbraid 7/0.2mm |
| Termination Type | Miniature or Standard Plug |
| Max. Temperature | 250°C |
| Min. Temperature | -50°C |
| Standards Met | ANSI |