



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## JIS Type K 1.9kg Pull Button Magnet Thermocouple, PFA Insulated Cable with Stainless Steel Over-Braid Terminating in Miniature Blue



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Specifications****Specifications**

<b>Product Code</b>	XE-0624-001
<b>General Description</b>	1.9kg Pull Button Magnet Thermocouple
<b>Sensor Type</b>	K JIS
<b>Magnet</b>	19.1mm Diameter / 12.7mm height - Max pull: 1.9 KG
<b>Cable Length</b>	1m
<b>Cable Type</b>	PFA Insulated Cable with stainless steel overbraid 7/0.2mm
<b>Termination Type</b>	Miniature Plug
<b>Max. Temperature</b>	250°C
<b>Min. Temperature</b>	-50°C
<b>Standards Met</b>	JIS (Japanese Industrial Standards)