

ANSI Crocodile Clip Thermocouple with Glassfibre Stainless Steel Overbraided Cable - Type K



[Click here to view product](#)

Product Code	XE-0628
General Description	Crocodile Clip Thermocouple with Glassfibre Stainless Steel Overbraided Cable
Thermocouple Type	K
Thermocouple Junction	Grounded measuring junction
Cable Length	1 or 3 metre
Cable Type	stainless steel over-braided PFA Teflon®
Cable Termination	Miniature Plug
Max. Temperature	250°C
Min. Temperature	-50°C
Standards Met	ANSI MC96.1