



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Fabricated Grounded Thermocouple in Stainless Steel Tube, Glassfibre stainless steel Overbraided cable - Type K JIS - 4.76mm diameter x 13mm length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-0649-001
Sensor Type	316 stainless steel tube in various sizes
Thermocouple Type	K JIS
Thermocouple Junction	Grounded at tip for fast response
Cable Length	2m
Cable Type	Glassfibre stainless steel overbraided cable
Cable Termination	Bare Wire Tails
Max. Temperature	350°C
Min. Temperature	-60°C
Standards Met	JIS (Japanese Industrial Standards)
Dimensions	4.76mm Diameter x 13mm Length