



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## **Fabricated Grounded Thermocouple in Stainless Steel Tube, Glassfibre stainless steel Overbraided cable - Type K JIS - 4.76mm diameter x 25mm length**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Specifications****Specifications**

<b>Product Code</b>	XE-0650-001
<b>Sensor Type</b>	316 stainless steel tube in various sizes
<b>Thermocouple Type</b>	K JIS
<b>Thermocouple Junction</b>	Grounded at tip for fast response
<b>Cable Length</b>	2m
<b>Cable Type</b>	Glassfibre stainless steel overbraided cable
<b>Cable Termination</b>	Bare Wire Tails
<b>Max. Temperature</b>	350°C
<b>Min. Temperature</b>	-60°C
<b>Standards Met</b>	JIS (Japanese Industrial Standards)
<b>Dimensions</b>	4.76mm Diameter x 25mm Length