

ANSI Fabricated Grounded Thermocouple in Stainless Steel Tube, Glassfibre stainless steel Overbraided cable - Type K



[Click here to view product](#)

Product Code	XE-0632
General Description	316 stainless steel tube in various sizes
Thermocouple Type	K
Thermocouple Junction	Grounded at tip for fast response
Cable Length	2m
Cable Type	Glassfibre stainless steel overbraided cable
Cable Termination	Bare Wire Tails
Max. Temperature	+350°C
Min. Temperature	-60°C
Standards Met	ANSI MC96.1